


SPECIFICATIONS

GeoPyc[®] 1365



REPRODUCIBILITY

When sample volume is at least 25% of sample holder volume: Typically $\pm 1.1\%$

SAMPLE CHAMBER

12.7-mm (0.50 in.) ID	Sample volumes between 0.3 cm ³ and 0.8 cm ³
19.1-mm (0.75 in.) ID	Sample volumes between 0.8 cm ³ and 2.4 cm ³
25.4-mm (1.00 in.) ID	Sample volumes between 2.4 cm ³ and 5.3 cm ³
38.1-mm (1.50 in.) ID	Sample volumes between 5.3 cm ³ and 13 cm ³
50.8-mm (2.00 in.) ID	Sample volumes between 13 cm ³ and 25 cm ³

PHYSICAL

Height	11 in.
Width	22 in.
Depth	15 in.
Weight	42 lbs

ELECTRICAL POWER

Electrical Power	95 VA
Voltage	85 to 265 VAC
Frequency	47 to 63 Hz

SPECIFICATIONS

GeoPyc[®] 1365

ENVIRONMENT

Temperature	15 °C to 35 °C
Humidity	20% to 80% non-condensing

**Due to continuous improvements, specifications are subject to change without notice.*