

BeVision D1 Specification:

Testing parameter	Material
Particle size distribution	Suspension, emulsion, dry powder
General	BeVision D1
Typical measurement time	Less than 10min
CCD	Imaging speed 120frames/s
Measuring speed	10000 pcs/min
Size	
Size range	Dry 30-10000um
Accuracy	≤1% (GBRM D50)
Repeatability	≤1% (GBRM D50)
Magnification	9-300 times
System	
	220VAC
Dimension	610mm x 300mm x 440mm (L x W x H)
Weight	26kg
Computer specification	
Computer interface	At least a USB2.0 port required
Operation system	Windows XP, Windows 7,8 or 10
Hardware specification	500GB HD, CPU: I7 2600 or above; Memory: 4G or above; Mainboard: With PCI-E interface, long and short card slot. Monitor resolution: 1680*1050 or above